

## CHAPTER 2

# Growth and the Environment of *Schizosaccharomyces pombe*

Janni Petersen<sup>1,3</sup> and Paul Russell<sup>2,3</sup>

<sup>1</sup>Flinders University, Flinders Centre for Innovation in Cancer, School of Medicine, FMST, Bedford Park, Adelaide, South Australia 5042, Australia; <sup>2</sup>Department of Cell and Molecular Biology, The Scripps Research Institute, La Jolla, California 92037

Here, we summarize the composition and uses of *Schizosaccharomyces pombe* media and discuss key issues for consideration in the generation of *S. pombe* cultures. We discuss the concept of “culture memory,” in which the growth state and stress experienced by a strain during storage, propagation, and starter culture preparation can alter experimental outcomes at later stages. We also describe the triggers that are widely used to manipulate signaling through the environment sensing pathways.

### INTRODUCTION

---

*Schizosaccharomyces pombe* cells are rod-shaped and grow by tip extension until they attain a critical cell length whereupon growth ceases and mitosis is initiated (Mitchison 1957). This highly reproducible cell size control is just one of the cellular growth and morphogenesis pathways that are regulated by environmental cues. The Sty1/Spc1 MAP kinase stress response and the target of rapamycin (TOR) pathways respond to changes in the extracellular environment to modulate functions as diverse as growth, division, the cytoskeleton, protein synthesis, transcription, and translation (Shiozaki and Russell 1995; Weisman 2004; Petersen 2009; Loewith and Hall 2011). The choice of media and culture conditions, therefore, heavily impacts the outcome of studies in virtually every arena of fission yeast biology.

### FISSION YEAST GROWTH MEDIA

---

Fission yeast has the ability to grow on a range of different media. These include complex rich medium, synthetic minimal medium, and different types of “sporulation” media used to promote sexual differentiation. Rich medium, known as yeast extract with supplements (YES) (see Recipes section below), is made with yeast extract and glucose with specific amino acids or nucleobases added as needed. YES supports fast growth rates (Table 1); hence, it is used for general nonselective growth. YES is also preferable for growth before long-term storage at  $-80^{\circ}\text{C}$  and for some experimental procedures, including protoplast generation. Synthetic media recipes are based on the original Edinburgh minimal media (EMM) developed by Murdoch Mitchison for his cell cycle analyses (Mitchison 1970). The totally defined nature of synthetic minimal media for logarithmic growth facilitates the reproduction of culture conditions between physiological experiments and laboratories, as well as the selection of a variety of auxotrophic markers on exogenous vectors or integrated DNA fragments.

---

<sup>3</sup>Correspondence: [janni.petersen@flinders.edu.au](mailto:janni.petersen@flinders.edu.au); [prussell@scripps.edu](mailto:prussell@scripps.edu)

Copyright © Cold Spring Harbor Laboratory Press; all rights reserved  
Cite this introduction as *Cold Spring Harb Protoc*; doi:10.1101/pdb.top079764

**TABLE 1.** Approximate generation times of healthy haploids according to medium (A) and nitrogen source (B), and cell dimensions according to ploidy (C)

A. Approximate generation times of healthy haploids		
Medium	Temperature	Generation time
YES	20°C	~8 h
	25°C	3 h
	28°C	2 h 45 min
	30°C	2 h 30 min
	32°C	2 h 10 min
	35°C	2 h
EMM2	20°C	~10 h
	25°C	4 h 30 min
	28°C	3 h 45 min
	30°C	3 h
	32°C	2 h 30 min
	35°C	2 h 20 min
B. Nitrogen-dependent generation times of healthy haploids		
Nitrogen source	Temperature	Generation time
EMM2 (ammonium)	28°C	3 h 45 min
EMMG (glutamic acid)	28°C	3 h 30 min
EMMP (proline)	28°C	4 h 30 min
MSL (arginine)	28°C	3 h 30 min
YES (complex)	28°C	2 h 45 min
C. Cell dimensions		
Ploidy	Cell length	Cell width (constant)
Haploid	12–14 $\mu\text{m}$	~3.5 $\mu\text{m}$
Diploid	20–24 $\mu\text{m}$	~4.5–5 $\mu\text{m}$

Changes to the nitrogen or carbon composition of synthetic minimal media are commonly used to manipulate growth and cell cycle controls and promote sexual differentiation.

### Carbon Sources in Growth Media

Although glucose is routinely used as the sole carbon source in growth medium, *S. pombe* can also use glycerol, sucrose, raffinose, and maltose (Schlanderer and Dellweg 1974). When fission yeast is grown on glucose alone, the levels of several respiratory enzymes are reduced (Poole and Lloyd 1973). Monitoring carbon dioxide production to gauge glycolytic flux has revealed that glycolysis accounts for most of the ATP production during aerobic growth in 3% glucose (Hamburger et al. 1977). In contrast, the catabolite repression that arises from propagation in glucose is lost when glycerol is the sole carbon source (Poole and Lloyd 1974), making growth in glucose highly reminiscent of the Warburg effect seen in cancer cells (Upadhyay et al. 2013). Removal of glucose while maintaining the provision of nitrogen can maintain cells in a  $G_0$ -like state for extended periods; this growth state is used to study *S. pombe* metabolic controls and ageing (Su et al. 1996).

### Nitrogen Sources in Growth Media

#### *Modified Versions of EMM*

Several modified versions of EMM can be selected according to the desired source of nitrogen (Table 2). Each distinct nitrogen source induces signature fluxes through the TOR and Sty1 growth control pathways to alter the architecture and flux of most cellular processes. The most widely used synthetic medium is EMM2 (see Recipes section below) (Nurse 1975), which employs the excellent nitrogen source ammonium chloride. The switch from the sodium acetate and sodium dihydrogen orthophosphate of the original EMM medium to 15 mM potassium hydrogen phthalate and disodium hydrogen

**TABLE 2.** Nitrogen sources in modified versions of EMM

Medium (EMM(X))	Nitrogen source
EMM2	Ammonium chloride
EMMG	L-Glutamic acid
EMMP	L-Proline
EMMSer	L-Serine
EMMUr	Uracil

orthophosphate was introduced by Paul Nurse to reduce clumping of *wee1-50* (also known as *cdc9-50*) mutant cells (Nurse 1975); however, as cells were generally healthier in EMM2, it became the medium of choice for most studies. EMM2 was further modified by Fantes and Nurse (1977) to assess the impact of nutrient variation on cell cycle controls, and different nitrogen sources were substituted for ammonium chloride. EMM2 without nitrogen is called EMM–N (see Recipes section below). (Note that EMM–N is commercially available from several providers, including Formedium in the U.K.) Replacement of the ammonium chloride of EMM2 with glutamic acid, proline, uracil, or serine generated the alternative media EMMG (see Recipes section below), EMMP (see Recipes section below), EMMUr (see Recipes section below), and EMMSer (see Recipes section below), respectively (Fantes and Nurse 1977). Of these media, EMMG, and EMMP have been the most widely used (Petersen and Nurse 2007). Because glutamic acid is a good nitrogen source and proline a very poor one, a switch between the two reveals the impact of changing the “quality” of the nitrogen source upon the cell cycle, stress, and TOR signaling networks (Petersen 2009). EMMP is also useful for cell cycle analysis, because cells are born in G<sub>1</sub> phase and initiate growth before committing to the cell cycle at START (Fantes and Nurse 1978). This contrasts with all other media, in which passage of START precedes cytokinesis.

Confusingly, the name EMMG was later used to describe an EMM2 variant in which the ammonium is replaced with 5.91 mM glutamate to support sporulation of diploids (Moreno et al. 1991). This led to the use of the name PMG (Pombe Minimal Glutamate) to describe the medium in which the ammonium of EMM2 is replaced with 20 mM glutamic acid (Sabatinos and Forsburg 2010). Throughout this collection, we return to the original nomenclature of EMMG for medium with glutamic acid. There is no mention in this collection of the medium with 5.91 mM glutamate described by Moreno et al. (1991). Wherever EMMG is used, it refers to medium with 20 mM glutamic acid.

EMMG is often a medium of choice when propagating certain auxotrophs. There appear to be two transporters for basic amino acids, of which one is inhibited by ammonium ions (Fantes and Creanor 1984). Growth on EMMG therefore has the advantage that both transporters can be used for uptake of basic amino acids. Empirical observations of a set of aromatic auxotrophic mutants show that EMMG supports much better growth than EMM2 in these mutants as well (Strauss 1979). For most researchers, however, EMMG is used because *ura4-d18* mutants are fitter when grown with EMMG than with EMM2, presumably because uracil uptake exploits the same channels. This negative impact of ammonium ions on the uptake of basic amino acids becomes a major issue when performing mass spectrometric analyses using SILAC (stable isotope labeling by amino acids in cell culture), in which strains are engineered to rely entirely upon the import of exogenous arginine and lysine that incorporate heavy carbon or nitrogen isotopes at key positions. Critically, the addition of lysine inhibits arginine uptake by the ammonium insensitive channel. As both lysine and arginine must be imported from the medium for SILAC to be successful, the high levels of ammonium in EMM2 either block growth altogether or generate an excessive lag phase before growth takes off. One solution is to reduce the ammonium levels in the EMM2 used for SILAC from 96 to 24 mM (Bicho et al. 2010); another is to switch to EMMG (Carpy et al. 2015).

### *Sporulation Media*

Minimal sporulating liquid or agar (MSL/MSA) medium (see Recipes section below) was developed by Egel et al. (1994) to facilitate the efficient induction of sexual differentiation in fission yeast. Like EMM-based media, it supports logarithmic growth; however, it employs arginine as its sole nitrogen

source. Egel selected arginine because a switch from MSL containing arginine to MSL without it, induced sexual differentiation far more rapidly and efficiently than a similar switch from EMM to EMM–N. Furthermore, when the background fluorescence of a fluorescent fusion protein obscures the specific signal in EMM2 medium, a switch from EMM2 to MSL often abolishes background fluorescence to reveal the “true” location of fusion-protein. A further reduction in background fluorescence in live cell imaging can be achieved by using a filter rather than autoclaving to sterilize medium. MSL and MSA are also useful for hygromycin and G418/kanamycin selection, which cannot be applied in EMM2. In contrast, sporulation agar (SPA) medium (see Recipes section below), in which no nitrogen is present, invokes an abrupt nitrogen starvation shock for cells (Leupold 1970). Cells plated on SPA either mate or arrest division, whereas they undergo one or two rounds of division when plated on MSA. MSA therefore produces efficient mating of most strains. However, spontaneous ascus breakdown is faster and more synchronous following mating on SPA than MSA, making SPA a medium of choice for mating strains for tetrad analysis.

### Supplementation of Minimal Media with Amino Acids and Nucleobases

For growth on solid media, 200 mg/L of the relevant supplements (e.g., leucine, uracil, histidine, adenine, lysine, arginine, or histidine) can be added. Note that amino acid supplements are often added to a final concentration of 225 or 250 mg/L; however, this high level of amino acids stems from Paul Nurse’s observation that survival of a mutant in *cdc2* required addition of 250 mg/L leucine. Consequently, 250 mg/L supplements are commonly used, although 150 mg/L is sufficient. It is notable that variance in size at division (i.e., variance in passage through the G<sub>2</sub> phase of the cell cycle) is lowest at 150 mg/L, rather than 250 mg/L (J Petersen and P Russell, unpubl.).

Amino acid stock solutions are often prepared as 7.5 mg/mL, however, both adenine and uracil have solubility issues. Adenine should be made as a 5-mg/mL stock solution and uracil solutions will need to be heated (e.g., by microwave) to dissolve the solid before adding to medium. The inconvenience and evaporation arising from repeated heating of the stock solution has prompted many to use the much more expensive uridine for ad hoc supplementation of media at the bench and only use uracil when making up large batches of medium.

#### *Mutations in ade6<sup>+</sup>*

Mutations in the *ade6<sup>+</sup>* gene are widely used in *S. pombe* genetics because a precursor in the adenine biosynthesis pathway, P-ribosylaminoimidazole, accumulates in the mutants and oxidizes to a red pigment (Reaume and Tatum 1949). Consequently, supplementation with a limiting level of adenine (10 µg/mL) supports growth yet reveals the red pigmentation (Schar and Kohli 1994). *ade6-M210* colonies appear dark red while *ade6-M216* colonies appear pink (Schar and Kohli 1994). Colonies harboring the *ade6-705* nonsense mutation are a deep crimson when not complemented by the *sup3-5* tRNA<sup>Ser</sup> mutation, but appear white when the opal suppressor allows reading through the stop mutation to produce a functional phosphor-ribosylaminoimidazole carboxylase enzyme (Hofer et al. 1979).

This color difference has been used for colony color/sectoring assays in genetic screens, for example, to detect the loss of a nonessential mini-chromosome. The suppression of *ade6<sup>+</sup>* mutations by mini-chromosome-borne sequences, which generates white colonies, is lost with the mini-chromosome. This produces neighboring red colonies or reveals a red sector within a colony when plates are supplemented with 10 µg/mL adenine (Niwa et al. 1986, 1989; Takahashi et al. 1994; Allshire et al. 1995).

### Phloxin B Staining

Phloxin B (2',4',5',7'-Tetrabromo-4,5,6,7-tetrachlorofluorescein disodium salt or cyanosine) is a red dye that is used as vital stain when added to YES; see yeast extract with supplements (YES) and phloxin B (see Recipes section below). Dead cells are unable to exclude phloxin B and thus appear red. This is especially useful for identifying colonies consisting of haploid or diploid cells, as the latter are more likely to die and therefore stain more darkly with phloxin B (Tange and Niwa 1995); thus, the overall

colony will be a darker pink. Phloxin B can also be useful for identifying colonies of temperature-sensitive lethal mutants because of increased cell death at the restrictive temperature.

## 5-FOA Selection

The orotidine 5'-phosphate decarboxylase encoded by the *ura4*<sup>+</sup> gene decarboxylates 5-fluoroorotic acid (5-FOA) to generate 5-fluorouracil, which blocks DNA replication. Therefore, strains from which the *ura4*<sup>+</sup> gene has been deleted or mutated are resistant to 5-FOA, making Yeast Extract with Supplements (YES) and 5-Fluoroorotic Acid (5-FOA) (see Recipes section below) an excellent reagent for positive selection for the loss of *ura4*<sup>+</sup> marker.

## STRAIN HUSBANDRY

---

### Storage

The best method for long-term storage of fission yeast is in 25% glycerol at  $-80^{\circ}\text{C}$ . After 2 d of growth in 10 mL of YES to generate a dense culture, cells are resuspended in 1 mL of YES containing 25% glycerol and transferred to a cryotube. If a plasmid that relies upon the selection of an auxotrophic marker must be maintained during storage, the strain should be grown with appropriate selection in minimal medium before resuspension in YES with 25% glycerol for storage at  $-80^{\circ}\text{C}$ . To reisolate strains from  $-80^{\circ}\text{C}$ , a sterile pipette tip is used to transfer a small lump of the frozen cell suspension to a “spot” on a plate (while taking care not to let the stock thaw). The cells are streaked out to generate single colonies whose genotypes should be checked before use.

For short-term storage, cells can be maintained at  $4^{\circ}\text{C}$  as patches on agar plates wrapped in plastic film for 2–3 wk. However, storage at  $4^{\circ}\text{C}$  can increase the frequency of diploidization and the emergence of phenotype-reversing mutations. It is therefore essential to reisolate the stored strain to a single colony and check the relevant genotypes (see below) each time it is used from a patch at  $4^{\circ}\text{C}$ .

### Checking Genotypes by Replica Plating

Replica plating is an efficient way to determine genotypes after storage. The plate containing the cells is inverted onto velvet or filter paper stretched over a replica plating block, and a low, even pressure is applied to ensure that a portion of each colony is retained on the paper after the plate is removed. The procedure is then repeated with dry blank plates to transfer cells from each colony on the paper to each blank plate. The addition or omission of appropriate supplements or drugs in these replica copies is used to reveal the genotype of each colony. Cells with a *ura4*<sup>-</sup> genotype divide several times after replica plating onto plates lacking uracil, making it necessary to sequentially replicate these colonies on two consecutive days to clearly differentiate uracil prototrophic from auxotrophic colonies.

## GROWTH STATE


---

Controlling cell physiology lies at the core of successful, reproducible experimentation with fission yeast. There are five principal impacts upon growth state: (1) changes in amino acid provision or requirement; (2) cell density; (3) the status of the starter culture used to inoculate the main culture; (4) the timing of inoculation; and (5) the level of aeration.

### Auxotrophic/Prototrophic Status and Amino Acid Provision

Auxotroph markers are invaluable tools for selecting for the maintenance of vectors or integration of foreign DNA sequences into the genome. When auxotrophs are propagated, supplements must be added to the medium to support cell growth. The addition of certain amino acids alone can radically change cell physiology (Weisman et al. 2005; Bonfils et al. 2012; Han et al. 2012; Kim et al. 2013; Xu





**FIGURE 1.** (A) Wild-type cells grown in EMMG with or without 150 mg/L leucine added. Cell lengths at division ( $\pm$ SEM, 200 cells) are shown. (B) Calcofluor staining of wild-type cells grown in medium (EMM2, EMMG, or EMMP) to a cell density of  $1.5 \times 10^6$  cells/mL. (C) Pheromone secretion by wild-type  $h^-$  cells. A diploid tester strain was mixed 100:1 with haploid  $h^-$  cells on MSA. Iodine staining was used to visualize the halo of sporulating diploid cells around the pheromone-producing haploid cells.

et al. 2013). For example, leucine addition transiently blocks mitotic commitment (Hartmuth and Petersen 2009) and alters cell length at division during steady-state growth (Fig. 1A). Therefore, it is essential to ensure that all strains in a given experiment have identical auxotrophic genotypes or, preferably, are all prototrophs. Amino acid dropout mixes that are commonly used for culturing *Saccharomyces cerevisiae* strains in synthetic growth media are not typically used for fission yeast, although they might be useful for isolating and culturing mutants that have caused unexpected auxotrophic phenotypes. If cell length at division is to be measured and the provision of amino acids cannot be avoided, supplementation with 150 mg/L of the relevant amino acid is recommended, as the standard deviation of cell length at division is lowest at this concentration of supplements (J Petersen and P Russell, unpubl.).

## Starvation-Induced Memory from Starters to Main Cultures

A widespread assumption, which may stem from work with budding yeast, is that because a new yeast cell is generated in one generation, a culture at any cell density will be the same irrespective of when it was inoculated. An extension of this misconception is the assumption that diluting a culture from a saturated stationary phase culture into the early log phase the evening before or the morning of an experiment will generate a culture from which meaningful conclusions can be drawn. However, this is not always the case. There is clear evidence of “culture memory” in *S. pombe*.

At the crudest level, the behavior of the spindle pole body shows that cells take at least two cell cycles to fully mature (Grallert et al. 2004). More generally, drop/spot test assays derived from cultures that were not in exponential phase can often produce misleading or irreproducible results. The ability to synchronize cultures by elutriation can be influenced by culture density. Inoculation from a dense culture the evening before elutriation often generates a size-selected population whose synchrony declines rapidly after the first division. Finally, the dynamic properties of microtubules vary according to culture age, time of inoculation and age of preculture (A Grallert, pers. comm.).

Changes in the heterochromatin that alter the global transcription landscape account for culture memory (Yamanaka et al. 2013). Alarming, once these modifications are established they may be stably inherited for many generations over several weeks of continual culture. There can be no stronger warning for taking extreme care to follow good practice in strain maintenance, propagation, and culture than this clear demonstration of culture memory.

## Physiological Experiments Using Liquid Cultures

The ability to grow *S. pombe* in defined minimal media has been a cornerstone of fission yeast research. Importantly, culture history (discussed above) and the density that a correctly inoculated culture achieves at the point of analysis are equally important to achieving reproducible results. For most studies, growth until early to mid-log phase is recommended:  $1 \times 10^6$ – $5 \times 10^6$  cells/mL in EMM2 and  $1 \times 10^6$ – $1 \times 10^7$  cells/mL in YES. Cell physiologies change significantly if cultures are grown beyond the upper ranges of these cell densities. For cell cycle analyses, comparisons should be made between cultures with identical biomass rather than cell densities. The optical density (OD) at 595 nm is a useful indirect measure of culture biomass (mass accumulation). Each spectrophotometer should be calibrated with cell number determined using a hemocytometer, but typically for wild-type haploid strains,  $OD_{595} = 0.1$  equates to  $\sim 1.5$ – $2 \times 10^6$  cells/mL.

## Culture Setup

The physiological state of a culture is influenced by the level of oxygenation; therefore, best practice dictates that a freshly isolated colony is inoculated into a 50-mL conical flask containing 5–10 mL of the same medium that will be used for the main culture. However, for particularly sick mutants, inoculating the starter culture in YES can be preferable. Following overnight shaking at the appropriate temperature, the starter culture should be diluted to early-exponential phase ( $\sim 1 \times 10^6$  cells/mL) if needed, and left to grow exponentially during the day. Alternatively, a freshly isolated colony can be inoculated into a nonshaking preculture which will then need to be grown with aeration and no starvation for at least 30 h before the experiment.

The main culture should be inoculated from the exponentially grown culture into a flask large enough to ensure the liquid is below the first line of graduation. Furthermore, when comparing two cultures, it is important to standardize aeration speed as well as culture and flask volume.

The following formula can be used to calculate the amount of starter culture needed to inoculate the final overnight culture:

$$\text{Volume (mL)} = \frac{(\text{Volume}_{\text{final}}[\text{mL}] \times \text{OD}_{\text{final}}) / 2^X}{\text{OD}_{\text{starter culture}}},$$

where  $X$  = number of generations = number of hours/generation time. If the starter culture is not in

exponential phase, one generation is subtracted to account for the time required for cells to reach optimum growth rate.

### Generation Time

The generation times for haploid wild-type *S. pombe* strains are listed in Table 1. Diploid strains have similar generation times. If the generation time of a mutant strain is unknown, it can be easily calculated from two cell number measurements (cells/mL) within exponential growth phase ( $1 \times 10^6$ – $5 \times 10^6$  cells/mL in minimal medium)

Generation time:  $T = \ln(2)/K$

$$K = \frac{\ln([\text{Cell number}_{T2}/\text{Cell number}_{T1}]^*)}{\text{Hours}(T2 - T1)}$$

\*For wild-type and mutant strains with normal coupling of growth and cell division, cells/mL at Time 2 (T2) and Time 1 (T1) can be replaced by OD<sub>595</sub> measurements at T2 and T1 as follows:

$$K = \frac{\ln(\text{OD}_{T2}/\text{OD}_{T1})}{\text{Hours}(T2 - T1)}$$

### Temperature

Wild-type fission yeast cells can grow between 18°C and 37°C, although growth at 37°C or below 20°C is far from ideal. However, this breadth of temperature range has been invaluablely exploited in the use of temperature-sensitive mutations that maintain function at 25°C but not at 35°C–37°C, or cold-sensitive mutations that maintain function above 30°C but not below 20°C.

### Nutrient Availability and Cell Physiology

Starvation arrests cell cycle progression to initiate a differentiation program that improves viability (Egel et al. 1990). When cultures reach stationary phase in most media, depletion of glucose will become the limiting factor and cells will arrest in G<sub>2</sub>. In contrast, cells arrest in G<sub>1</sub> when specifically starved of nitrogen (a response that prepares them for sexual differentiation). When the nitrogen source is changed, proliferation continues but at a rate that is specific for each nitrogen source (Table 1; Fig. 1B; Fantès and Nurse 1977). This distinction in the physiological response to growth on different nitrogen sources makes it difficult to compare experiments performed in different media (Table 1). For example, the *plp1.S402E* mutation has no impact upon cell length at division (a measure of the duration of G<sub>2</sub> phase) in EMM2, but decreases it by 1.1 μm in EMMG (Petersen and Hagan 2005; Petersen and Nurse 2007).

Changes in the flux through the “stress response” Sty1/Spc1 MAP kinase cascade are revealed by changes in activating phosphorylation of the Sty1/Spc1 by the MAPK kinase Wis1, or in the nuclear import of Sty1-GFP (Shiozaki and Russell 1995; Gaits et al. 1998; Toone and Jones 2004). Importantly, the steady-state “nonstressed” Sty1/Spc1 and TOR activities differ depending on the medium of choice (Shiozaki and Russell 1995; Petersen and Nurse 2007). A final point that is often ignored is that the slight differences in cell density arising from the gradual depletion of nutrient supply during growth alter cell physiology. Consequently, cell length at division of wild-type prototrophs increases as cell density/biomass increases (Fig. 1A). It is therefore essential to compare length at division at identical biomass (OD measurement) rather than cell number.

### Sources of Unwanted Stress in Liquid Cultures

Most studies inadvertently induce stress by harvesting cells through centrifugation. Exposure to the 800 g routinely used to harvest cells in a microfuge enhances Sty1/Spc1 signaling (Shiozaki et al. 1998), making cell isolation by filtration the only option for some experiments. However, it is important to



ensure that filters are not overloaded or dried out during filtration, and that the conical flask and medium receiving the filtered cells are prewarmed to the appropriate temperature.

## GROWTH RATE MEASUREMENTS

---

### Monitoring Optical Density

Monitoring OD at 595 nm provides a measurement of growth rate or biomass accumulation. A parallel estimate of cell numbers over time can be used to determine whether the exponential increase in cell growth and cell division displays the linear correlation of wild-type strains under steady state growth conditions. Importantly, in wild-type cells this coupling between cell growth and division is altered by environmental changes including starvation. Cell counting with a hemocytometer is a simple way to monitor cell number increases, although more accurate measurements can be attained using a Coulter counter.

### Assessing Fitness with Spot/Drop Tests

Testing the ability of 10-fold serial dilutions of cultures to form colonies on agar plates provides a powerful assessment of gene function and the impact of multiple genetic interactions upon cell viability. It is particularly useful for detecting subtle differences in growth rate or viability. It is essential that the cultures used to inoculate such spot/drop tests have been grown under identical physiological conditions (in the same way and to the same cell density), and are diluted 10-fold without any prior centrifugation or filtration. While spot/drop tests have traditionally been end point analyses that catalog colony formation after 2–3 d, technologies developed for genome-wide analyses can be used to monitor colony formation by time-lapse microscopy (Banks et al. 2012). This new approach has proven highly informative, indicating that checking plates on a daily basis reveals early differences that are lost when colonies saturate toward the end of growth.

## NITROGEN STARVATION

---

Nitrogen starvation can be used to arrest the cell cycle in  $G_1$  (Egel 1971; Nurse 1975); however, it is important to remember that cell metabolism is altered by this starvation. If both mating types are present, nitrogen starvation will induce sexual differentiation. If sexual differentiation is to be monitored in liquid cultures, increasing aeration through the use of baffled conical flasks will considerably enhance the efficiency of differentiation. (The expression of agglutinins ensures cells aggregate for mating in the most turbulent of conditions.)

### Nitrogen Starvation to Enhance Phenotype Penetrance in Deletion Spores

Heterozygous diploids, in which one copy of a gene of interest has been replaced with a selectable marker, can be sporulated by switching a log-phase culture into nitrogen-free medium. Preferential growth of the deletant spores over the nondeletants can be ensured when the spore mix is inoculated into liquid cultures by selecting for or against the gene encoded in the deletion cassette. Although such spores will have low levels of the protein encoded by the deleted gene, some protein and mRNAs may be inherited from the parental zygote to confound phenotype analysis. This residual protein can be eliminated by germination in nitrogen-free EMM2 (EMM2–N). The lack of nitrogen prompts autophagy to scavenge nitrogen for protein synthesis by turning over existing proteins, including any residual parental protein. Thus, the germinating spore depletes the protein of interest so that the subsequent addition of nitrogen provides true insight into the “null” phenotype.

## Starvation-Induced Pseudohyphal Growth

When nitrogen deprivation is combined with a good carbon source at a remote location, pseudohyphal growth produces elaborately branched pseudohyphae that invade deep into solid medium to scavenge for nutrients (Amoah-Buahin et al. 2005).

## HALO ASSAYS

---

Mating pheromones are secreted by fission yeast upon nitrogen starvation. Egel et al. (1994) developed a semiquantitative plate assay for monitoring pheromone secretion. A diploid strain that only sporulates upon exposure to exogenous pheromone is used to monitor the gradient of pheromone produced by the strain of interest. Iodine staining reveals a halo of dark asci around individual cells (Fig. 1C). The diameter of the halo is indicative of the level of pheromone secreted. In principle, halo assays—in which drug-soaked filter papers are placed on lawns of cells—should be as effective in fission yeast as they are in budding yeast for the assessment of cell growth in response to small molecules. If a small molecule blocks growth, a clearing zone of no growth will be observed (Gassner et al. 2007).

## GROWTH IN RESPONSE TO STRESS AGENTS

---

Although the Sty1/Spc1 MAPK pathway is activated in response to stress (Toone and Jones 2004), it is rarely appreciated that a low level activation of the Sty1/Spc1 pathway under steady-state, or “non-stressed,” conditions makes a significant contribution to the timing of cell division (Shiozaki and Russell 1995). Thus, it may be more appropriate to consider the Sty1/Spc1 pathway as a homeostasis, rather than stress, pathway. Cellular stresses simply enhance the signaling through the pathway, which is already setting a level of response for a particular steady state. Stimulation of the Sty1/Spc1 pathway by a wide range of agents promotes a remarkably sophisticated set of changes to help the cell deal with each individual insult with a specific response.

### Osmotic Stress

Either 1.2 M sorbitol or 0.65–1.0 M KCl is added to plates to invoke an osmotic stress response; however, it is important to note that, unlike sorbitol, KCl also induces salt ion stress. For osmotic stress in liquid cultures, 0.6–0.9 M KCl results in progressively increased activation of Spc1/Sty1.

### Temperature

Survival after brief exposure to higher temperatures (42°C) is commonly used to assess tolerance to heat stress. Furthermore, steady-state growth at 37°C is stressful, and mutants in the Sty1/Spc1 signaling pathway cannot grow at this temperature.

### Nutrients

The reduction in the quality of nutrient supply in a shift from EMMG to EMMP stimulates Sty1/Spc1 signaling (Fantes and Nurse 1977; Hartmuth and Petersen 2009). The generation time in medium using proline as a nitrogen source is longer compared with medium exploiting glutamic acid (Table 1), and cells are born in the G<sub>1</sub> rather than G<sub>2</sub> phase of the cell cycle.

### Heavy Metals

Exposure to heavy metals, including copper and cadmium, stimulates Sty1/Spc1 signaling to promote survival. CdCl<sub>2</sub> or CdSO<sub>4</sub> is added to medium at a final concentration of 0.2 mM (Ortiz et al. 1992),

and  $\text{CuCl}_2$  or  $\text{CuSO}_4$  to a concentration of 16 nM (Labbe et al. 1997). Metalloid arsenite (0.5 M) also activates Spc1/Sty1 (Rodriguez-Gabriel and Russell 2005).

## Oxidative Stress

Oxidative stress can be generated with several different chemicals. The most commonly used oxidative stress agent is hydrogen peroxide ( $\text{H}_2\text{O}_2$ ). Hydrogen peroxide is highly unstable; therefore, it is best to freeze small aliquots from a new bottle and discard these ~1 wk after thawing. The concentration of hydrogen peroxide can be modified to suit the experimental goals, but as a general guideline a final concentration of 0.2 mM  $\text{H}_2\text{O}_2$  is sufficient to cause rapid and robust activation of Sty1/Spc1. The growth of wild-type cells is inhibited by 0.5–1.0 mM  $\text{H}_2\text{O}_2$  in plates containing YES and lower concentrations in plates containing EMM2. *Tert*-butyl hydroperoxide, also known as *t*-BuOOH or under the trade name Luperox, can also be used at concentrations of ~1 mM to cause oxidative stress.

## DNA Damage

Fission yeast has been an excellent model organism for studying cellular responses to DNA damage and replication inhibition. A large variety of DNA damaging agents have been used for these studies. Ionizing radiation (IR) is often used to cause double-strand breaks (DSBs), which are among the most lethal forms of DNA damage. The exact protocols vary depending on the specific type of  $\gamma$  irradiator. For laboratories that do not have access to a  $\gamma$  irradiator, the glycopeptide antibiotic bleomycin (~5 mU/mL) and the related compound phleomycin (~10  $\mu\text{g}/\text{mL}$ ) are effective substitutes for IR. The precise mechanism of action for these genotoxins is unknown, but it is believed to involve a generation of free radicals in close proximity to DNA. The quinolone alkaloid camptothecin (CPT), which is a DNA topoisomerase I (topo I) inhibitor, can also be used to create DSBs. In this case, DSBs arise when replication forks encounter cleavage complexes formed by DNA, CPT, and topo I, thereby converting single-strand DNA breaks into one-ended DSBs. CPT is typically used at concentrations between 0.1 and 10  $\mu\text{M}$ . The DNA alkylating agent methyl methanesulfonate (MMS) is used in agar plates at concentrations between 0.001% and 0.01%. Toxicity from methylated DNA primarily arises from replication fork stalling and collapse. Ultraviolet (UV) radiation causes several types of DNA damage, most notably cyclobutane–pyrimidine dimers (CPDs) and 6-4 photoproducts (6-4PPs). UV exposure can be performed with a UV cross-linker (254 nm), originally designed for fixing purified nucleic acids to membranes. In this event, plates must be dry and cells present in a monolayer. Hydroxyurea (HU), an inhibitor of ribonucleotide reductase, causes replication forks to stall by depleting deoxyribonucleotides. HU is used at concentrations in the range of 0.1–10 mM and upon extended exposure can lead to DNA strand breaks.

## GROWTH CONDITIONS FOR EFFICIENT MATING

---

Malt extract agar (MEA) medium (see Recipes section below) and EMM–N can be used to invoke sexual differentiation, but SPA and MSL both support more efficient and more rapid mating (Gutz and Doe 1973). Zygotes form after 2 d on EMM2–N but overnight on SPA. When setting up crosses, there are a number of points to bear in mind. It is important to remember that *S. pombe* can recycle the nitrogen of added amino acids to use as a nitrogen source for the synthesis of all amino acids. Thus, when mating auxotrophs, it is important to supply only the amino acids required to support growth of the two parental strains (i.e., do not simply use a mixture of the five most commonly required amino acids) and to supply them at a minimal level (5–40  $\mu\text{g}/\text{mL}$ ). The nitrogen starvation response and, consequently, mating efficiency declines with every microgram of amino acid supplied. Although wild-type cells can mate from colonies or patches on plates stored at 4°C, it is important to remember that all discussions around the benefits of working with healthy cells for proliferating cultures apply to

the optimization of mating efficiency. Cells entering the stationary phase invoked by nutrient exhaustion on old plates mate poorly, while those entering from actively growing, fresh, healthy colonies are much more efficient. When a mating persistently fails it can be useful to change to another mating medium. If this is still unsuccessful, then mixing equal numbers of cells from mid-log phase liquid cultures and washing them in nitrogen-free medium before plating them on the surface of the mating plate can persuade the most reluctant of partners to form zygotes. A final point to remember is that mating is blocked above 33°C and probably most efficient at 30°C.

## RECIPES

---

### *EMM2*

Reagent	Quantity (for 1 L)	Final concentration
Ammonium chloride	5 g	93.5 mM
Potassium hydrogen phthalate	3 g	14.7 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15.5 mM
Glucose	20 g	111 mM
Salt stock (50×)	20 mL	1×
Vitamins (1000×)	1 mL	1×
Minerals (10,000×)	0.1 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare 1 L of EMM2 by combining the reagents above. To prepare solid medium, include agar (20 g/L). Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). (For liquid medium, low pressure and a short cycle is essential to avoid caramelization of glucose and breakdown of vitamins and minerals. For live cell imaging, filter-sterilized medium has lower background fluorescence than autoclaved medium. Use a 0.22- $\mu$ m pore size for filtration.) Store at 4°C. If desired, add 12  $\mu$ M thiamine (364  $\mu$ L of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

### *EMMG*

Reagent	Quantity (for 1 L)	Final concentration
L-Glutamic acid, monosodium salt	3.38 g	20 mM
Potassium hydrogen phthalate	3 g	14.7 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15.5 mM
Glucose	20 g	111 mM
Salt stock (50×)	20 mL	1×
Vitamins (1000×)	1 mL	1×
Minerals (10,000×)	0.1 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare 1 L of EMMG by combining the reagents above. To prepare solid medium, include agar (20 g/L). Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). (For liquid medium, low pressure and a short cycle is essential to avoid caramelization of glucose and breakdown of vitamins and minerals. For live cell imaging, filter-sterilized medium has lower background fluorescence than autoclaved medium. Use a 0.22- $\mu$ m pore size for filtration.) Store at 4°C. If desired, add 12  $\mu$ M thiamine (364  $\mu$ L of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

### EMM-N

Reagent	Quantity (for 1 L)	Final concentration
Potassium hydrogen phthalate	3 g	14.7 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15.5 mM
Glucose	20 g	111 mM
Salt stock (50×)	20 mL	1×
Vitamins (1000×)	1 mL	1×
Minerals (10,000×)	0.1 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare 1 L of EMM-N (EMM2 without nitrogen) (Nurse 1975) by combining the reagents above. To prepare solid medium, include agar (20 g/L). Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). Store at 4°C. If desired, add 12 μM thiamine (364 μL of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

### EMMP

Reagent	Quantity (for 1 L)	Final concentration
L-Proline	2.3 g	20 mM
Potassium hydrogen phthalate	3 g	14.7 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15.5 mM
Glucose	20 g	111 mM
Salt stock (50×)	20 mL	1×
Vitamins (1000×)	1 mL	1×
Minerals (10,000×)	0.1 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare 1 L of EMMP by combining the reagents above. To prepare solid medium, include agar (20 g/L). Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). (For liquid medium, low pressure and a short cycle is essential to avoid caramelization of glucose and breakdown of vitamins and minerals. For live cell imaging, filter-sterilized medium has lower background fluorescence than autoclaved medium. Use a 0.22-μm pore size for filtration.) Store at 4°C. If desired, add 12 μM thiamine (364 μL of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

### EMMSer

Reagent	Quantity (for 1 L)	Final concentration
L-Serine	2.1 g	20 mM
Potassium hydrogen phthalate	3 g	14.7 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15.5 mM
Glucose	20 g	111 mM
Salt stock (50×)	20 mL	1×
Vitamins (1000×)	1 mL	1×
Minerals (10,000×)	0.1 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare 1 L of EMMSer by combining the reagents above. To prepare solid medium, include agar (20 g/L). Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). (For liquid medium, low pressure and a short cycle is essential to avoid caramelization of glucose and breakdown of vitamins and minerals. For live cell imaging, filter-sterilized medium has lower background fluorescence than autoclaved medium. Use a 0.22-μm pore size for filtration.) Store at 4°C. If desired, add 12 μM thiamine (364 μL of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).



**EMMUr**

Reagent	Quantity (for 1 L)	Final concentration
Uracil	1.12 g	10 mM
Potassium hydrogen phthalate	3 g	14.7 mM
Na <sub>2</sub> HPO <sub>4</sub>	2.2 g	15.5 mM
Glucose	20 g	111 mM
Salt stock (50×)	20 mL	1×
Vitamins (1000×)	1 mL	1×
Minerals (10,000×)	0.1 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare 1 L of EMMUr by combining the reagents above. To prepare solid medium, include agar (20 g/L). Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). (For liquid medium, low pressure and a short cycle is essential to avoid caramelization of glucose and breakdown of vitamins and minerals. For live cell imaging, filter-sterilized medium has lower background fluorescence than autoclaved medium. Use a 0.22- $\mu$ m pore size for filtration.) Store at 4°C. If desired, add 12  $\mu$ M thiamine (364  $\mu$ L of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

**Malt Extract Agar (MEA) Medium**

For 1 L of MEA (Gutz and Doe 1973), mix 30 g of malt extract with 20 g of agar and bring to 1 L with H<sub>2</sub>O. Sterilize by autoclaving at 10 psi for 15 min. Store at 4°C.

**Minerals (10,000×)**

Reagent	Quantity (for 1 L)	Final concentration
Boric acid	5 g	80.9 mM
MnSO <sub>4</sub>	4 g	23.7 mM
ZnSO <sub>4</sub> · 7H <sub>2</sub> O	4 g	13.9 mM
FeCl <sub>2</sub> · 6H <sub>2</sub> O	2 g	7.40 mM
Molybdic acid	0.4 g	2.47 mM
KI	1 g	6.02 mM
CuSO <sub>4</sub> · 5H <sub>2</sub> O	0.4 g	1.60 mM
Citric acid	10 g	47.6 mM

Store at 4°C.

**Minimal Sporulating Liquid or Agar (MSL/MSA) Medium**

Reagent	Quantity (for 1 L)	Final concentration
Glucose	10 g	55.5 mM
Arginine hydrochloride	2 g	9.5 mM
KH <sub>2</sub> PO <sub>4</sub>	1 g	7.3 mM
NaCl	0.1 g	1.7 mM
MgSO <sub>4</sub> · 7H <sub>2</sub> O	0.2 g	0.81 mM
CaCl <sub>2</sub> · 6H <sub>2</sub> O	0.15 g	0.68 mM
Minerals (10,000×)	0.1 mL	1×
Vitamins (1000×)	2 mL	1×
H <sub>2</sub> O	to 1 L	

Prepare minimal sporulation liquid (MSL) medium (Egel et al. 1994) by combining the reagents in the list above. To prepare solid minimal sporulating agar (MSA) medium, include 20 g/L of agar. Sterilize by autoclaving at 10 psi for 10 min (for liquid medium) or 15 min (for solid medium). (For liquid medium, low pressure and a short cycle is essential to avoid caramelization of glucose and breakdown of vitamins and minerals. For live cell imaging, filter-sterilized medium has lower

background fluorescence than autoclaved medium. Use a 0.22- $\mu$ m pore size for filtration.) Store at 4°C. If desired, add 12  $\mu$ M thiamine (364  $\mu$ L of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

**Salt Stock (50 $\times$ )**

Reagent	Quantity (for 1 L)	Final concentration
MgCl <sub>2</sub> · 6H <sub>2</sub> O	52.5 g	0.26 M
CaCl <sub>2</sub> · 2H <sub>2</sub> O	0.735 g	4.99 mM
KCl	50 g	0.67 M
Na <sub>2</sub> SO <sub>4</sub>	2 g	14.1 mM

Store at 4°C.

**Sporulation Agar (SPA) Medium**

Reagent	Quantity (for 1 L)	Final concentration
Glucose	10 g	55.5 mM
KH <sub>2</sub> PO <sub>4</sub>	1 g	7.3 mM
Vitamins (1000 $\times$ )	1 mL	1 $\times$
Agar	20 g	20 g/L
H <sub>2</sub> O	to 1 L	

To prepare SPA medium (Munz and Leupold 1970), combine the ingredients listed above. Sterilize by autoclaving at 10 psi for 15 min. Store at 4°C. If desired, add 12  $\mu$ M thiamine (364  $\mu$ L of a filter-sterilized stock solution of 10 mg/mL on H<sub>2</sub>O) after autoclaving to fully repress expression from *nmt1*-derived promoters (Maundrell 1990).

**Vitamins (1000 $\times$ )**

Reagent	Quantity (for 1 L)	Final concentration
Pantothenic acid	1 g	4.20 mM
Nicotinic acid	10 g	81.2 mM
Inositol	10 g	55.5 mM
Biotin	10 mg	40.8 $\mu$ M

Store at 4°C.

**Yeast Extract with Supplements (YES)**

Reagent	Amount to add (1 L)
Glucose	30 g
Yeast extract	5 g
Adenine	0.2 g
Uracil	0.2 g
Histidine	0.2 g
Leucine	0.2 g
Lysine	0.2 g
Agar (for solid medium only)	20 g
H <sub>2</sub> O	to 1 L

Prepare 1 L of YES by combining the reagents listed above. Sterilize by autoclaving at 10 psi for 15 min. Store at 4°C. (Note that the yeast extract contains sufficient thiamine to repress *nmt1* promoters.)

*Yeast Extract with Supplements (YES) and 5-Fluoroorotic Acid (5-FOA)*

Reagent	Amount to add (1 L)
Glucose	60 g
Yeast extract	10 g
Adenine	0.4 g
Uracil	0.1 g
Histidine	0.4 g
Leucine	0.4 g
Lysine	0.4 g
Agar (for solid medium only)	40 g
H <sub>2</sub> O	to 1 L

Prepare 1 L of 2× YES by combining the reagents listed above. Sterilize by autoclaving at 10 psi for 15 min. Dissolve 2 mg of 5-FOA (Sigma-Aldrich F5013) in 1 mL of sterile milliQ H<sub>2</sub>O. Mix 1:1 with 55°C 2× YES using a sterile stir bar. (5-FOA from Sigma-Aldrich is very expensive; cheaper reagents from other suppliers do not work as well. Calculate the quantities that you will require to make up the precise number of plates you need for your selection.) Store at 4°C. (Note that the yeast extract contains sufficient thiamine to repress *nmt1* promoters.)

*Yeast Extract with Supplements (YES) and Phloxin B*

Reagent	Amount to add (1 L)
Glucose	30 g
Yeast extract	5 g
Adenine	0.2 g
Uracil	0.2 g
Histidine	0.2 g
Leucine	0.2 g
Lysine	0.2 g
Agar (for solid medium only)	20 g
H <sub>2</sub> O	to 1 L

Prepare 1 L of YES by combining the reagents listed above. Sterilize by autoclaving at 10 psi for 15 min. Add 4 mL of Phloxin B stock solution (5 g/L, filter-sterilized) to 1 L of YES after autoclaving. Store at 4°C. (Note that the yeast extract contains sufficient thiamine to repress *nmt1* promoters.)

---

## ACKNOWLEDGMENTS

We thank K. Kowalczyk for technical assistance with the Halo assay. J.P. is supported by Manchester and Flinders University and a Cancer Research UK Senior Fellowship [C10888/A11178] and P.R. is supported by grants from the National Institutes of Health [GM059447, CA077325, CA117638, and P42ES010337].

---

## REFERENCES

- Allshire RC, Nimmo ER, Ekwall K, Javerzat JP, Cranston G. 1995. Mutations derepressing silent centromeric domains in fission yeast disrupt chromosome segregation. *Genes Dev* **9**: 218–233.
- Amoah-Buahin E, Bone N, Armstrong J. 2005. Hyphal growth in the fission yeast *Schizosaccharomyces pombe*. *Eukaryot Cell* **4**: 1287–1297.
- Banks AP, Lawless C, Lydall DA. 2012. A quantitative fitness analysis workflow. *J Vis Exp* doi: 10.3791/4018.
- Bicho CC, de Lima Alves F, Chen ZA, Rappsilber J, Sawin KE. 2010. A genetic engineering solution to the “arginine conversion problem” in stable isotope labeling by amino acids in cell culture (SILAC). *Mol Cell Proteomics* **9**: 1567–1577.
- Bonfils G, Jaquenoud M, Bontron S, Ostrowicz C, Ungermann C, De Virgilio C. 2012. Leucyl-tRNA synthetase controls TORC1 via the EGO complex. *Mol Cell* **46**: 105–110.

- Carpay A, Patel A, Tay YD, Hagan IM, Macek B. 2015. Nic1 inactivation enables stable isotope labeling with  $^{13}\text{C}_6$   $^{15}\text{N}_4$ -arginine in *Schizosaccharomyces pombe*. *Mol Cell Proteomics* 14: 243–250.
- Egel R. 1971. Physiological aspects of conjugation in fission yeast. *Planta* 98: 89–96.
- Egel R, Nielsen O, Weilguny D. 1990. Sexual differentiation in fission yeast. *Trends Genet* 6: 369–373.
- Egel R, Willer M, Kjaerulf S, Davey J, Nielsen O. 1994. Assessment of pheromone production and response in fission yeast by a halo test of induced sporulation. *Yeast* 10: 1347–1354.
- Fantes P, Nurse P. 1977. Control of cell size at division in fission yeast by growth-modulated size control over nuclear division. *Exp Cell Res* 107: 377–386.
- Fantes PA, Nurse P. 1978. Control of the timing of cell division in fission yeast. Cell size mutants reveal a second control pathway. *Exp Cell Res* 115: 317–329.
- Fantes PA, Creanor J. 1984. Canavanine resistance and the mechanism of arginine uptake in the fission yeast *Schizosaccharomyces pombe*. *J Gen Microbiol* 130: 3265–3273.
- Gaits F, Degols G, Shiozaki K, Russell P. 1998. Phosphorylation and association with the transcription factor Atf1 regulate localization of Spc1/Sty1 stress-activated kinase in fission yeast. *Genes Dev* 12: 1464–1473.
- Gassner NC, Tumble CM, Bock JE, Cotton N, White KN, Tenney K, St Onge RP, Proctor MJ, Giaever G, Nislow C, et al. 2007. Accelerating the discovery of biologically active small molecules using a high-throughput yeast halo assay. *J Nat Prod* 70: 383–390.
- Grallert A, Krapp A, Bagley S, Simanis V, Hagan IM. 2004. Recruitment of NIMA kinase shows that maturation of the *S. pombe* spindle-pole body occurs over consecutive cell cycles and reveals a role for NIMA in modulating SIN activity. *Genes Dev* 18: 1007–1021.
- Gutz H, Doe FJ. 1973. Two different  $h^+$  mating types in *Schizosaccharomyces pombe*. *Genetics* 74: 563–569.
- Hamburger K, Kramhohf B, Nissen SB, Zeuthen E. 1977. Linear increase in glycolytic activity through the cell cycle of *Schizosaccharomyces pombe*. *J Cell Sci* 24: 69–79.
- Han JM, Jeong SJ, Park MC, Kim G, Kwon NH, Kim HK, Ha SH, Ryu SH, Kim S. 2012. Leucyl-tRNA synthetase is an intracellular leucine sensor for the mTORC1-signaling pathway. *Cell* 149: 410–424.
- Hartmuth S, Petersen J. 2009. Fission yeast Tor1 functions as part of TORC1 to control mitotic entry through the stress MAPK pathway following nutrient stress. *J Cell Sci* 122: 1737–1746.
- Hofer F, Hollenstein H, Janner F, Minet M, Thuriaux P, Leupold U. 1979. The genetic fine structure of nonsense suppressors in *Schizosaccharomyces pombe*. I. *sup3* and *sup9*. *Curr Genet* 1: 45–61.
- Kim J, Song G, Wu G, Gao H, Johnson GA, Bazer FW. 2013. Arginine, leucine, and glutamine stimulate proliferation of porcine trophoblast cells through the MTOR-RPS6K-RPS6-EIF4EBP1 signal transduction pathway. *Biol Reprod* 88: 113.
- Labbe S, Zhu Z, Thiele DJ. 1997. Copper-specific transcriptional repression of yeast genes encoding critical components in the copper transport pathway. *J Biol Chem* 272: 15951–15958.
- Leupold U. 1970. Genetical methods for *Schizosaccharomyces pombe*. *Methods Cell Physiol* 4: 169–177.
- Loewith R, Hall MN. 2011. Target of rapamycin (TOR) in nutrient signaling and growth control. *Genetics* 189: 1177–1201.
- Maundrell K. 1990. Nmt1 of fission yeast. A highly transcribed gene completely repressed by thiamine. *J Biol Chem* 265: 10857–10864.
- Mitchison JM. 1957. The growth of single cells. I. *Schizosaccharomyces pombe*. *Exp Cell Res* 13: 244–262.
- Mitchison JM. 1970. Physiological and cytological methods for *Schizosaccharomyces pombe*. *Methods Cell Physiol* 4: 131–165.
- Moreno S, Klar A, Nurse P. 1991. Molecular genetic analysis of fission yeast *Schizosaccharomyces pombe*. *Methods Enzymol* 194: 795–823.
- Munz P, Leupold U. 1970. Characterization of ICR-170-induced mutations in *Schizosaccharomyces pombe*. *Mutat Res* 9: 199–212.
- Niwa O, Matsumoto T, Yanagida M. 1986. Construction of a mini-chromosome by deletion and its mitotic and meiotic behavior in fission yeast. *Mol Gen Genet* 203: 397–405.
- Niwa O, Matsumoto T, Chikashige Y, Yanagida M. 1989. Characterization of *Schizosaccharomyces pombe* minichromosome deletion derivatives and a functional allocation of their centromere. *EMBO J* 8: 3045–3052.
- Nurse P. 1975. Genetic control of cell size at cell division in yeast. *Nature* 256: 547–551.
- Ortiz DF, Kreppel L, Speiser DM, Scheel G, McDonald G, Ow DW. 1992. Heavy metal tolerance in the fission yeast requires an ATP-binding cassette-type vacuolar membrane transporter. *EMBO J* 11: 3491–3499.
- Petersen J. 2009. TOR signalling regulates mitotic commitment through stress-activated MAPK and Polo kinase in response to nutrient stress. *Biochem Soc Trans* 37: 273–277.
- Petersen J, Hagan IM. 2005. Polo kinase links the stress pathway to cell cycle control and tip growth in fission yeast. *Nature* 435: 507–512.
- Petersen J, Nurse P. 2007. TOR signalling regulates mitotic commitment through the stress MAP kinase pathway and the Polo and Cdc2 kinases. *Nat Cell Biol* 9: 1263–1272.
- Poole RK, Lloyd D. 1973. Oscillations of enzyme activities during the cell cycle of a glucose-repressed fission-yeast *Schizosaccharomyces pombe* 972h-. *Biochem J* 136(1): 195–207.
- Poole RK, Lloyd D. 1974. Changes in respiratory activities during the cell cycle of the fission yeast *Schizosaccharomyces pombe* 972h- growing in the presence of glycerol. *Biochem J* 144(1): 141–148.
- Reaume SE, Tatum EL. 1949. Spontaneous and nitrogen mustard-induced nutritional deficiencies in *Saccharomyces cerevisiae*. *Arch Biochem* 22: 331–338.
- Rodriguez-Gabriel MA, Russell P. 2005. Distinct signaling pathways respond to arsenite and reactive oxygen species in *Schizosaccharomyces pombe*. *Eukaryot Cell* 4: 1396–1402.
- Sabatino SA, Forsburg SL. 2010. Molecular genetics of *Schizosaccharomyces pombe*. *Methods Enzymol* 470: 759–795.
- Schar P, Kohli J. 1994. Preferential strand transfer and hybrid DNA formation at the recombination hotspot *ade6-M26* of *Schizosaccharomyces pombe*. *EMBO J* 13: 5212–5219.
- Schlenderer G, Dellweg H. 1974. Cyclid AMP and catabolite repression in yeasts, in *Schizosaccharomyces pombe* glucose lowers both intracellular adenosine 3': 5'-monophosphate levels and the activity of catabolite-sensitive enzymes. *Eur J Biochem* 49(1): 305–316.
- Shiozaki K, Russell P. 1995. Cell-cycle control linked to extracellular environment by MAP kinase pathway in fission yeast. *Nature* 378: 739–743.
- Shiozaki K, Shiozaki M, Russell P. 1998. Heat stress activates fission yeast Spc1/Sty1 MAPK by a MEKK-independent mechanism. *Mol Biol Cell* 9: 1339–1349.
- Strauss A. 1979. Isolation and genetical classification of aromatic amino acid auxotrophic mutants in *Schizosaccharomyces pombe*. *J Gen Microbiol* 113: 173–176.
- Su SS, Tanaka Y, Samejima I, Tanaka K, Yanagida M. 1996. A nitrogen starvation-induced dormant  $G_0$  state in fission yeast: The establishment from uncommitted  $G_1$  state and its delay for return to proliferation. *J Cell Sci* 109: 1347–1357.
- Takahashi K, Yamada H, Yanagida M. 1994. Fission yeast minichromosome loss mutants *mis* cause lethal aneuploidy and replication abnormality. *Mol Biol Cell* 5: 1145–1158.
- Tange Y, Niwa O. 1995. A selection system for diploid and against haploid cells in *Schizosaccharomyces pombe*. *Mol Gen Genet* 248: 644–648.
- Toone WM, Jones N. 2004. Stress responses in *S. pombe*. In *The molecular biology of Schizosaccharomyces pombe* (ed. Egel R), pp. 57–72. Springer, Berlin, Heidelberg.
- Upadhyay M, Samal J, Kandpal M, Singh OV, Vivekanandan P. 2013. The Warburg effect: Insights from the past decade. *Pharmacol Ther* 137: 318–330.
- Weisman R. 2004. The fission yeast TOR proteins and the rapamycin response: An unexpected tale. *Curr Top Microbiol Immunol* 279: 85–95.
- Weisman R, Roitburg I, Nahari T, Kupiec M. 2005. Regulation of leucine uptake by *tor1<sup>+</sup>* in *Schizosaccharomyces pombe* is sensitive to rapamycin. *Genetics* 169: 539–550.
- Xu B, Lee KK, Zhang L, Gerton JL. 2013. Stimulation of mTORC1 with L-leucine rescues defects associated with Roberts syndrome. *PLoS Genet* 9: e1003857.
- Yamanaka S, Mehta S, Reyes-Turcu FE, Zhuang F, Fuchs RT, Rong Y, Robb GB, Grewal SI. 2013. RNAi triggered by specialized machinery silences developmental genes and retrotransposons. *Nature* 493: 557–560.